


da Vinci Radical Cystectomy maintains the oncologic principles of open radical cystectomy via a minimally invasive approach.

Patient Positioning & Preparation


- Pad, secure and prep patient while in supine position (legs abducted, thighs at table-level)
- DVT prophylaxis
- Shave patient from costal margins to pubic bone
- Sterilely prepare and drape abdomen, penis, scrotum, perineum and peri-anal region
- Insert Foley catheter and rectal bougie
- Place patient in steep Trendelenburg position



Steep Trendelenburg Position (>20°)

Port Placement (Assumes a left-side assistant)


- **da Vinci Camera Port, 12 mm (Blue):** Place near umbilicus, 20 cm above symphysis pubis.
- **Right da Vinci Instrument Arm ①, 8 mm (Yellow):** Place 8 cm from the endoscope port along the line from the umbilicus to the right anterior superior iliac spine (ASIS).
- **Left da Vinci Instrument Arm ②, 8 mm (Green):** Place 8 cm from the endoscope port along the line from the umbilicus to the left ASIS.
- **3rd da Vinci Instrument Arm ③, 8 mm (Red):** Place 8 cm directly lateral to the right instrument port.
- **Assistant Port (A1), 12 mm (White):** Place two-finger breadths medial and cephalad to the left ASIS, and at least 8 cm from the left instrument port.
- **Assistant Port (A2), 12 mm (White):** Place 10 cm cephalad to the right instrument port.
- **Assistant Port (A3), 5 mm (Black):** Place 12 cm cephalad to the left instrument port.



Port Placement

MALE: Radical Cystectomy (Cystoprostatectomy)

Surgery to remove the bladder, prostate, seminal vesicles and lymph nodes




The green shading indicates anatomy to be removed.

FEMALE: Radical Cystectomy (Anterior Pelvic Exenteration)

Surgery to remove the bladder, urethra, uterus, cervix, ovaries, anterior wall of the vagina* and lymph nodes

*Spare vagina when appropriate.



Patient Cart Positioning

Place patient in steep Trendelenburg position (>20°) before rolling in patient cart. Align camera arm, camera arm set-up joint and center column with target anatomy (bladder). Position and dock patient cart directly between legs.

Docking

• Dock the camera arm first

Use port and arm clutch maneuvers to ensure camera is aligned with target anatomy. If 3rd instrument arm will be used, shift the camera arm set-up joint towards the patient's side with just one instrument arm.

• Dock the remaining instrument arms

Maximize spacing between all instrument arms.



4-Arm da Vinci System
(set-up assumes left-side assistant)

da Vinci Radical Cystectomy - Procedure Steps, Instrumentation & Accessories

Surgical Step	Recommended EndoWrist® Instrument	Camera/Scope & Ancillary Supplies
Isolate Ureters - Dissect the ureters distally to the level of the bladder; ligate and transect. Female: Identify and divide ovarian pedicles prior to ureteral isolation.	Monopolar Instrument <i>Hot Shears™</i> (Monopolar Curved Scissors) or Permanent Cautery Hook	0° scope; large <i>Hem-o-lok®</i> clips
Posterior Dissection - Incise the peritoneum between the bladder and rectum; mobilize the bladder and posterior prostate off the rectum and dissect to Denonvilliers' fascia. Female: Dissect between the uterus/vagina and bladder. May enter and excise anterior vaginal wall en bloc with bladder or spare vagina when appropriate.		30° scope (up) if needed
Lateral Dissection - Incise the peritoneum just lateral to the medial umbilical ligaments bilaterally. Develop the space between the bladder and lateral pelvic side wall using blunt dissection. Female: Complete division of the cardinal and uterosacral ligament bilaterally.		
Secure Bladder Pedicles - Ligate and transect the pedicles using an endoscopic linear stapler/cutter, or bipolar cautery and clips.		stapler/cutter or large <i>Hem-o-lok®</i> clips
Anterior Dissection - Drop the bladder (as is done in <i>da Vinci®</i> Prostatectomy) by incising the peritoneum anteriorly.		
Control of Prostatic Apex and Urethra - Incise the endopelvic fascia and develop the prostatic apex. Ligate the DVC and transect the anterior aspect of urethra.	Bipolar Instrument Maryland Bipolar Forceps or <i>PK™</i> Dissecting Forceps	DVC stitch: 0 Vicryl™ on CT-2 (cut to 6")
Prostatic Dissection (Male Only) - Use clips and/or bipolar cautery to secure pedicles. Nerve-Sparing (Select Patients): Use clips and cold scissors to secure pedicles.		large <i>Hem-o-lok®</i> clips
Hysterectomy and Bi-Lateral Salpingo-Oophorectomy (Female Only) - Retract the uterus superiorly, excise the uterus and remove the fallopian tubes, ovaries and cervix with the bladder specimen.	3rd Arm Options Fenestrated Bipolar Forceps <i>ProGrasp™</i> Forceps Cadiere Forceps Cobra Grasper	
Bag Specimen - Place the specimen in a retrieval bag and bring into upper quadrants of the abdomen. Irrigate and ensure hemostasis.		specimen bag
Reconstruct Vagina (Female Only) - Reconstruct vagina when needed.		0-Vicryl™ on CT-2
Bilateral Pelvic Lymphadenectomy - Lymphadenectomy generally includes: obturator nodes, external and common iliac nodes. Para-aortic dissection is possible with the <i>da Vinci® S™</i> System.		
Tag the Ureters/Pre-Place Anastomotic Stitch (Option 1) - Tag ureters for ready localization during open diversion. In cases of neobladder, pre-place posterior anastomotic stitches in urethra if performing the urethra-neobladder anastomosis extracorporeally.		3-0 Vicryl™ on SH, full-length
Urinary Diversion (Neobladder or Ileal Conduit) - Extract specimen, tunnel ureters and harvest bowel through the 6-8 cm periumbilical midline incision. Perform the urinary diversion and ureteroenteric anastomosis extracorporeally.	Other Instruments <i>Hem-o-lok®</i> Clip Applier, Large	
Urethra-Neobladder Anastomosis (Option 2) - Utilize pre-placed urethral anastomotic stitches if performing the urethra-neobladder anastomosis extracorporeally. If anastomosis is to be performed robotically, close the periumbilical incision, dock the patient cart and complete the anastomosis.		